

Site AED Inspections

#	AED Location	Pad Exp Dates	July				August				September								
			1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk					
1																			
2																			
3																			
4																			
5																			
6																			
7																			
8																			
9																			
10																			
#	AED Location	Pad Exp Dates	October				November				December								
1			1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk					
2																			
3																			
4																			
5																			
6																			
7																			
8																			
9																			
10																			

Suggested signaturred weekly walk by inspectors, checking for ready light or OK symbol or obvious alarms or problem indicator. All employees should look for and report any potential problems or malfunctions of their AEDs to their immediate supervisor. Visible, directional signs should in place; follow manufacturer's recommended AED testing guidelines. Any changes to AED units, location or usage, please contact Collier County EMS - 239-252-3740.

Site AED Inspections

#	AED Location	Pad Exp Dates	January				February				March						
			1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk			
1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	

#	AED Location	Pad Exp Dates	April				May				June						
			1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk	1st wk	2nd wk	3rd wk	4th wk			
1																	
2																	
3																	
4																	
5																	
6																	
7																	
8																	
9																	
10																	

Suggested signed weekly walk by inspections, checking for ready light or OK symbol or obvious alarms or problem indicator. All employees should look for and report any potential problems or malfunctions of their AEDs to their immediate supervisor. Visible, directional signs should be in place, follow manufacturers recommended AED testing guidelines. Any changes to AED units, location or usage, please contact Collier County EMS - 239-252-3740.